

# Enterella

## SX-ELLA Stent Colorectal



Easy deployment

Excellent flexibility

Atraumatic stent ends

Smooth alignment around flexures



## Basic description

### SX-ELLA Stent Colorectal – Enterella

Stent system for fluoroscopic insertion  
Self-expandable nitinol stent preloaded in a pull-back delivery system

## Indications

Temporary or definitive palliative treatment of malignant strictures affecting the rectum, sigmoid colon or distal part of descending colon  
Treatment of coloenteral and/or colovesical fistulas (covered stent)

## Contraindications

Bowel perforation  
Peritoneal abscess  
Severe ascites  
Severe coagulopathy  
Enteral ischemia

## Features / Benefits


- ✓ Atraumatic stent ends / Effective prevention of trauma to the bowel wall
- ✓ Excellent flexibility / Precise alignment around anatomical curves
- ✓ Radiopaque markers at both stent ends and at the midpoint / Improved visibility; accurate stent positioning
- ✓ Plastic repositioning loop at aboral stent end / Possibility of stent repositioning within five days after implantation
- ✓ Optional stent covering / Prevention of tumor in-growth; occlusion of fistulas

## Additional description

15 F (5.08 mm) delivery system available in a standard length of 95 cm (135 cm delivery system is available on request, covered stent requires 18 F (6.0 mm) delivery system).

Radiopaque, atraumatic olive tip of the delivery system.  
Recommended guide wire: 0.035" (0.89 mm) ultra stiff.


## Excellent flexibility



## Plastic repositioning loop



## Atraumatic stent ends




## Easy deployment



## Available sizes

REF No.	Nominal (unconstrained) diameter [mm]	Nominal (unconstrained) length [mm]
019-05-22-082	22	82
019-05-22-090	22	90
019-05-22-113	22	113
019-05-22-135	22	135
019-05-25-082	25	82
019-05-25-090	25	90
019-05-25-113	25	113
019-05-25-135	25	135
019-05-30-075	30	75
019-05-30-088	30	88
019-05-30-112	30	112
019-05-30-123	30	123
019-05-30-136	30	136



Distributor:

Manufacturer:

**ELLA-CS, s.r.o.**

Milady Horakove 504, 500 06 Hradec Kralove  
Czech Republic

Tel.: +420 495 279 111 Fax: +420 495 265 655  
E-mail: info@ellacs.eu URL: www.ellacs.eu