



Enterella

SX-ELLA Stent Pyloroduodenal



Enterella




Easy deployment

Excellent flexibility

Atraumatic stent ends

Smooth alignment around flexures



Basic description

SX-ELLA Stent Pyloroduodenal – Enterella

Stent system for fluoroscopic and / or endoscopic insertion
Self-expandable nitinol stent preloaded in a pull-back delivery system

Indications

Temporary or definitive palliative treatment of malignant strictures affecting gastric outlet and / or duodenum

Contraindications

Enteral ischemia
Bowel perforation
Peritoneal abscess
Severe ascites
Severe coagulopathy


Features / Benefits

- ✓ Atraumatic stent ends / Effective prevention of trauma to the gastric / duodenal wall
- ✓ Excellent flexibility / Precise conformation to the anatomical curves in a patient
- ✓ Radiopaque markers at both stent ends and at the midpoint / Improved visibility; accurate stent positioning
- ✓ Plastic repositioning loop at proximal stent end / Possibility of stent repositioning within five days following implantation
- ✓ Braided kink-resistant delivery system / Superior system introduction and stent deployment

Additional description

10.5 F (3.5 mm) delivery system available in a standard length of 210 cm allowing passage through the working channel of an endoscope
Radiopaque, atraumatic olive tip of the delivery system
Recommended guide wire: 0.035" (0.89 mm)


Kink-resistant delivery system




Yellow visual marker and radiopaque ring




Radiopaque markers




Plastic repositioning loop



Atraumatic stent ends



Easy deployment



Available sizes

REF No.	Nominal (unconstrained) diameter [mm]	Nominal (unconstrained) length [mm]
019-04-20-082	20	82
019-04-20-090	20	90
019-04-20-113	20	113
019-04-20-135	20	135
019-04-22-082	22	82
019-04-22-090	22	90
019-04-22-113	22	113
019-04-22-135	22	135
019-04-25-082	25	82
019-04-25-090	25	90
019-04-25-113	25	113
019-04-25-135	25	135

Distributor:

Manufacturer:

ELLA-CS, s.r.o.

Milady Horakove 504, 500 06 Hradec Kralove
Czech Republic

Tel.: +420 495 279 111 Fax: +420 495 265 655
E-mail: info@ellacs.eu URL: www.ellacs.eu